

PRODUCT SPECIFICATION

Product: SALVEQUICK DETECTABLE
ADHESIVE BANDAGE (PLASTER)

Reference code: 6754CAP



1. General description

Detectable, sterile adhesive bandages of two sizes, medium and fingertip, packed in a cardboard refill. The sterile adhesive bandages are covered by a paper cap. The bandages are detectable by most common metal detectors (ferro/non-ferro). The refill fits into the Cederroth (Salvequick) dispenser for adhesive bandages.

The adhesive bandage consists of a backing, evenly spread with an adhesive coating. An absorbent pad is fixed centrally to the dressing leaving a margin of adhesive surface around the pad. The adhesive is covered with a release-coated paper /foil that consists of two parts partly overlapping each other. When taking the adhesive bandage out from the package, the lower part of the cover will remain behind. The cap is still on the dressing. One half of the dressing can now be applied on the skin using one hand only. After removing the cap and the upper part of the cover the adhesive bandage is ready for use.

2. Intended use

The sterile adhesive bandages are intended to be used as a mechanical barrier and for absorption of exudates from minor wounds.

3. Technical description

Product dimensions

Medium	Adhesive bandage:	72x19 mm
	Wound pad:	24x12 mm
Fingertip	Adhesive bandage:	66x39 mm
	Wound pad:	28x12 mm

Material

Backing	Blue, soft Polyurethan (PU) foil
Adhesive	Transparent polyacrylate (hypoallergenic)
Woundpad	Laminate of dry laid, needle punched non-woven with a non-stick wound contact layer and a metal detectable aluminum foil.
Facing (cover)	Siliconised paper and PE-foil
Wrapping	Plain coldseal paper + paper cap

4. Sterilization

Sterilization method: Gamma irradiation

5. Composition

Refill

<u>Ref #</u>	<u>Product/Material</u>	<u>Pcs</u>
45167121	Medium Plastic adhesive bandage with cap	15
45167140	Fingertip Plastic adhesive bandage with cap	15
Dimensions (l x w x h) mm:		30 x 87 x 100

Consumer package

Number of refills in consumer package:	6
Dimensions (l x w x h) mm:	100 x 150 x 90

Transport package

Number of consumer packages in transport package:	10
Dimensions (l x w x h) mm:	520 x 160 x 190
Weight (kg):	2,2

6. Labeling

Product (refill)

Cederroth address, telephone- and faxnumber
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Symbols: CE +NB#, STERILE, R, DO NOT REUSE, LOT, USE BY

Transport package

Cederroth address, telephone- and faxnumber, web-address
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Quantity

Consumer package

Cederroth address, telephone- and fax-number, web- and mail-address
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Quantity
 Symbols: CE +NB#, LOT, USE BY

7. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

8. Accessories

The following accessories can be used with the product;

Ref #: 1907 Salvequick Plaster Dispenser



Authorized by: Teija Ålander
Quality Assurance Manager

Date: 2008-10-30